JOINTLY
BUILDING BRIDGES

C+TBA
austria

Allograft C+TBA

- DE GEBRAUCHS- UND FACHINFORMATION
- EN PACKAGE LEAFLET: INFORMATION FOR THE USER
- FR MODE D'EMPLOI ET INFORMATIONS TECHNIQUES
- ES PROSPECTOY INFORMACIÓN TÉCNICA DEL PRODUCTO
- NL GEBRUIKS- EN VAKINFORMATIE
- DK INDLÆGGSSEDEL OG PRODUKTRESUMÉ
- NO PAKNINGSVEDLEGG OG FAGINFORMASJON
- FI PAKKAUSSELOSTE JA VALMISTEYHTEENVETO
- IT FOGLIETTO ILLUSTRATIVO E RIASSUNTO DELLE CARATTERISTICHE DEL PRODOTTO
- GR ΕΛ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΚΑΙ ΕΙΔΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ
- SE BIPACKSEDEL OCH PRODUKTRESUMÉ

DE

EN

FR

ES

NL

DK

NO

FI

IT

GR

SE

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Allograft C⁺TBA

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EN

DE

PACKAGE LEAFLET: INFORMATION FOR THE USER

Please read carefully

EN

1. Name of the tissue product

Allograft C+TBA

Allogenic bone tissue of human origin, lyophilised; C+TBA

FR

2. Composition

Cancellous portion of human bone tissue

Cortical portions of human bone tissue

Corticospungious portions of human bone tissue

ES

3. Description

Bone graft of human origin, lyophilised, in cancellous, cortical or cortico-cancellous form as blocks, diaphyses, granules, rings, plates, cubes and cylinder.

NL

4. Clinical particulars

4.1 Indications for Use

Allograft C+TBA is suitable for use on toddlers, children, adolescents and adults for transplantation as a replacement for lost and insufficient tissue and for filling up and stabilising bone defects in orthopaedics, trauma surgery as well as oral and maxillofacial surgery. Use in further specialist disciplines is conceivable.

DK

NO

FI

IT

4.2 Dosage, method and duration of administration

Aseptic conditions are essential when opening the Allograft C+TBA packaging and for subsequent handling. If the implant is contaminated during the surgical procedure, it must be discarded.

Before application in patients, molded bodies in particular should be rehydrated in a suitable physiological medium for at least 10 minutes (e.g. physiological saline).

Rehydration with the following quantities of reconstitution solution is recommended for the product Granules + NaHyA. After thorough mixing, a moldable texture is formed, which can be adjusted to the desired consistency by adding the reconstitution solution drop by drop (~0.1 ml). If the reconstitution quantity is under or exceeded, the described consistency may not form or may be lost. Similarly, excessive liquid at the augmentation site, e.g. in the case of heavy bleeding in sinus floor or socket augmentation, can impair the handling of the bone grafting material during application and

GR

SE

DE

its adhesion to the defect. It is recommended to remove excess fluid at the defect site before application.

EN

FR

ES

NL

DK

NO

FI

IT

GR

SE

Particle size	Volume specification of tissue product	Volume of reconstitution solution
< 0,25 mm	0,25 ml	0,15 ± 0,05 ml
	0,5 ml	0,3 ± 0,1 ml
	1 ml	0,7 ± 0,1 ml
	2 ml	1,4 ± 0,2 ml
	5 ml	3,5 ± 0,2 ml
	10 ml	7,0 ± 0,2 ml
0,25-1 mm	0,25 ml	0,2 ± 0,1 ml
	0,5 ml	0,4 ± 0,1 ml
	1 ml	0,8 ± 0,1 ml
	2 ml	1,6 ± 0,2 ml
	5 ml	4 ± 0,2 ml
	10 ml	8 ± 0,2 ml
1-2 mm	0,25 ml	0,3 ± 0,1 ml
	0,5 ml	0,6 ± 0,1 ml
	1 ml	1,2 ± 0,2 ml
	2 ml	2,5 ± 0,2 ml
	5 ml	5,0 ± 0,2 ml
	10 ml	10 ± 0,2 ml

The area to be filled must always be completely filled with bone substitute material.

The final wound closure must be complete and saliva-proof. Wound closure should be free of tension. The allogenic bone tissue must be stabilised in order to avoid any micro-movement.

The material is dense to X-rays, and therefore radiological control presents no problem.

Allograft C-TBA is intended for single patient use only. Unused product residues must be discarded after the completion of the surgical procedure.

Preparation of recipient (patient):

1. First the surgeon has to remove necrotic or inflamed tissue. Any foreign materials and sequestra that are present also have to be removed and the host bed has to be cleaned.
2. The bone surface has to be freshened and has to bleed slightly.

Advice for the patient:

If the graft is not completely structurally supported, always warn the patient that he/she is not allowed to put excessive strain on the respective body site through carrying heavy weights over a longer period.

4.3 Contraindications

Necrosis in the surgical area is a contraindication for Allograft C+TBA. Following the complete removal of infected bone tissue the use of Allograft C+TBA should be considered carefully due to limited healing rate. Persons with known hypersensitivity and/or allergy to the ingredient sodium hyaluronate must be excluded from treatment with Granules + NaHyA. For patients with malignant tumor diseases, a contraindication exists due to the suspicion that NaHyA exhibits a cell-proliferative property at certain cell receptors.

4.4 Special warnings and precautions for use

Allograft C+TBA remains sterile as long as the packaging is not opened or damaged. If the packaging is externally damaged, it must be assumed that the product is non-sterile.

When implanting multiple products, the patient's body weight should be taken into account. In adults weighing at least 60 kg, a maximum of 20 ml of tissue (equivalent to, for example, 20 units for products with 1 ml each) should be used per surgery. For lower body weights, a proportionally reduced amount of products should be used accordingly. The application in host bed with poor blood circulation is subject to closer examination and consideration by the attending clinician due to poor healing rates.

Although no case of immunological rejection was observed in the longstanding use of human bone allografts, the possibility of such a rejection cannot be ruled out completely. After the initial testing of each donor, the batch processing takes place in accordance with all technical safety specifications, whereby various donors can be combined.

4.5 Interactions with other medicinal products and other forms of interactions

None known

4.6 Pregnancy and lactation

There is no adequate data on use during pregnancy and lactation. The risks that are generally linked to a surgical procedure including drug administration, blood loss etc., must be especially taken into careful consideration for pregnant and lactating patients. Each case should be considered carefully.

DE

EN

FR

ES

NL

DK

NO

FI

IT

GR

SE

DE

4.7 Effects on ability to drive and use machines:

Allograft C+TBA does not affect the ability to drive or use machines. During the surgical procedure, medicinal products such as analgesics may be used which affect the ability to drive or use machines. The package leaflet and instructions for use of the corresponding medicinal products must be taken into account.

EN

FR

ES

NL

DK

NO

FI

IT

GR

SE

4.8 Undesirable effects

Frequencies are defined using the following conventions: Very common ($\geq 1/10$); common ($\geq 1/100$, $< 1/10$); uncommon ($\geq 1/1,000$, $< 1/100$); rare ($\geq 1/10,000$, $< 1/1,000$); very rare ($< 1/10,000$); unknown (frequency cannot be estimated from the available data)

As with all bone grafts, the use of surgery makes it difficult to assess whether side effects, such as infections, complications during healing, or loosening of the graft are caused by the medicinal product or by the surgery.

A retrospective analysis demonstrated the safety of human lyophilised allografts. No microbial or viral infections were caused by the allograft.

Side effects known to occur due to the surgery include haematomas and pain.

Donors of bone material that is processed into Allograft C+TBA are subject to strict selection criteria with respect to donor selection and serological testing. A validated inactivation procedure is used for the processing. Despite these extensive safety measures, it is not possible to exclude with absolute certainty the transmission of infectious diseases due to the presence of unknown or undetectable pathogens.

4.9 Overdosage

Not applicable

5. Pharmacological properties

5.1 Pharmacodynamic properties

The pharmacodynamic effect of the medicinal product is due to its osteoconductive property.

A summary of clinical studies with human allogeneic bone tissue showed its good incorporation, a lack of rejection and good remodelling.

A retrospective analysis of the efficacy and safety of human lyophilised allografts demonstrated their safe use and osteo-conductivity.

Radiation has no negative effect on stability, healing success or incorporation of the transplant.

5.2 Pharmacokinetic properties

Not applicable

5.3 Preclinical safety data

As this medicinal product has been used for years in man, no preclinical animal studies have been performed.

6. Pharmaceutical Particulars

6.1 List of inactive ingredients

None

6.2 Incompatibilities

Incompatibilities were not identified during clinical studies. Concomitant medication with the diclofenac substance group may impair bone growth.

6.3 Shelf life

in the closed packaging – see imprint on the label and/or the packaging. Do not use the medicinal product after expiry of its shelf life.

6.4 Special precautions for storage

Store at 5-30 °C.

Store in the original packaging, in order to protect the contents from humidity. Allograft C-TBA is transported at an ambient temperature and protected from moisture.

6.5 Nature and contents of container

The Allograft C-TBA products have double packaging. The primary packaging (granules: vial closed with a lid or blister with sealing film, molded bodies: sealed peel-off packaging) keeps the product sterile and is suitable for use in the surgery. The secondary packaging is a sealed peel-off packaging.

6.6 Special precautions for disposal and other handling instructions

The content of the package is intended for single-patient use only. Any remaining residues should be discarded. Disposal must be carried out according to local regulations.

The product must not be re-sterilised.

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Only available by prescription

Transplant certificate:

A transplant certificate is attached to each product or can be found at the following link <https://www.ctba.at/begleitschein>.

Immediate postoperative adverse drug reactions (especially infections) and their accompanying circumstances must be listed on the attached notification sheet. The transplant certificate is used to record and evaluate clinical results and side effects and should be carefully filled out and returned after completion of the treatment, but no later than 6 months post-op. Particular emphasis is also placed on the transmission of verbal information, which should reflect experiences in the clinical use of the products or contain suggestions for their further development.

General Advice:

Allograft C-TBA can be ordered by licensed doctors or pharmacists.

Shipping usually happens via courier.

In case of any serious adverse reactions observed during and after clinical application which may be linked to Allograft C-TBA notify tissue establishment Cells+Tissuebank Austria gemeinnützige GmbH

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