

REBSTOCK®
Precision made in Germany

NEW!



STERILE CRANIO FIXATIONSSETS
STERILE CRANIAL FIXATION SETS

/// www.rebstock.de

Artikel-Nr.
Cat. No.



S50-31-004

Cranio Fixationsset, steril, Standard, 2-Loch
Cranial Fixation Sterile Set, standard, 2-holes

bestehend aus:
consisting of:

- 50-31-102 3x** Platte, gerade, 2-Loch mit 12 mm Steg
plate, straight, 2-holes with 12 mm bar
- 50-31-352 6x** Schraube, Kreuzschlitz, 1,6 mm x 4 mm, selbstbohrend
screw, cross head, 1.6 mm x 4 mm, self drilling



S50-31-005

Cranio Fixationsset, steril, Standard, 4-Loch
Cranial Fixation Sterile Set, standard, 4-holes

bestehend aus:
consisting of:

- 50-31-104 3x** Platte, gerade, 4-Loch mit 6 mm Steg
plate, straight, 4-holes with 6 mm bar
- 50-31-352 10x** Schraube, Kreuzschlitz, 1,6 mm x 4 mm, selbstbohrend
screw, cross head, 1.6 mm x 4 mm, self drilling



S50-31-006

Cranio Fixationsset, steril, Bohrlochplatte
Cranial Fixation Sterile Set, Burr Hole Cover

bestehend aus:
consisting of:

- 50-31-139 1x** Bohrlochplatte, ø 15 mm
burr hole plate, ø 15 mm
- 50-31-102 2x** Platte, gerade, 2-Loch mit 12 mm Steg
plate, straight, 2-holes with 12 mm bar
- 50-31-352 10x** Schraube, Kreuzschlitz, 1,6 mm x 4 mm, selbstbohrend
screw, cross head, 1.6 mm x 4 mm, self drilling

REBSTOCK INSTRUMENTS GmbH

In Weiheräcker 7
78589 Dürbheim-Tuttlingen
GERMANY

Fon +49 7424 / 98230-30
Fax +49 7424 / 98230-50
info@rebstock.de

/// www.rebstock.de