

NEONATAL / PEDIATRIC

*Neonatal & Pediatric Surgery
Specialised Newborn Baby's Surgery*





Designed specifically to meet

*the needs of Neonates - when it comes
to infant and baby surgical techniques.*

NEONATAL & PEDIATRIC SURGERY

Index

Diethrich DeBakey Bulldog	4
Diethrich Non-Traumatic Bulldog	5
Micro Bulldog, Weldon (Glover) Mini Bulldog	6
DeBakey Neonatal Ultralight Clamp	7
Cooley Neonatal Ultralight Clamp	9
Castaneda Neonatal Clamp	10
DeBakey Ring Handle Bulldog, Mini-Gregory Clamp	11
Contractor Bailey Baby, Thorax Retractor	12
Selman Cooley Pediatric Clamp	15
Cooley Pediatric Vascular Clamp, Cooley Anastomosis Clamp	16
Cooley Pediatric Clamp	17
Cooley Spoon-Shaped Clamp, Cooley Derra Clamp	20
Cooley Clamp	21
Kartchner Carotid Clamp, DeBakey Clamp	22
Cooley Kitzmiller Pediatric Clamp, Cooley Martinez Clamp	23
DeBakey Pediatric Clamp	24
DeBakey Vascular Clamp, DeBakey Glover Coarctation Clamp	28
DeBakey Peripheral Vascular Clamp, Eriksson-Stille DeBakey Carotid Clamp	29
Gregory Profunda Clamp, DeBakey Bainbridge Clamp	30
Leland-Jones Clamp, Castaneda Pediatric Clamp	31
Instrumentation Supplies	32

Legend

All graphics are original size (100%) if not otherwise noted.

Grey catalogue numbers are only available on request.

Blue marked articles (RT) are also available in titanium.

Overall length depends on jaw and handle style.

All measurements, pressures and angles are typical and may vary due to production.



CASTANEDA jaw/tip
1 row of microscopic teeth



COOLEY jaw/tip
2 rows of microscopic teeth



DEBAKEY jaw/tip
1x2 rows of microscopic teeth


This Reuchlen catalog is effective with October 2019.

The content is exclusively for information purposes, can be changed without prior notice at any time and can not be seen as our obligation. The August Reuchlen GmbH grants no warranty or guarantee with regard to correctness and accuracy of the content in this catalog.

Instrument design can be subject to change or modification without notice.

© Copyright August Reuchlen GmbH®, Tuttlingen, Germany


DIETRICH DEBAKEY BULLDOG

RC 2045

standard pressure ca. 80g	RC 2043 53 / 12mm	RC 2049 56 / 15mm	RC 2045 63 / 23mm
increased pressure ca. 120g	RC 2043-120 53 / 12mm	RC 2049-120 56 / 15mm	RC 2045-120 63 / 23mm

Diethrich DeBakey
Cross-Action Bulldog





RC 2047

standard pressure ca. 80g	RC 2044 52 / 10mm	RC 2048 54 / 14mm	RC 2047 57 / 20mm
increased pressure ca. 120g	RC 2044-120 52 / 10mm	RC 2048-120 54 / 14mm	RC 2047-120 57 / 20mm


Diethrich DeBakey
Cross-Action Bulldog


DIETRICH NON-TRAUMATIC BULLDOG

 RC 2123	standard pressure ca. 80g	RC 2123 52 / 8mm	RC 2125 56 / 12mm	RC 2127 60 / 16mm	RC 2129 64 / 20mm
	reduced pressure ca. 60g	RC 2123-060 52 / 8mm	RC 2125-060 56 / 12mm		
	increased pressure ca. 120g	RC 2123-120 52 / 8mm	RC 2125-120 56 / 12mm		RC 2129-120 64 / 20mm


Diethrich Non-Traumatic
Cross-Action Bulldog



 RC 2124	standard pressure ca. 80g	RC 2124 52 / 8mm	RC 2126 54 / 12mm	RC 2128 56 / 16mm	RC 2130 60 / 20mm
	reduced pressure ca. 60g	RC 2124-060 52 / 8mm	RC 2126-060 54 / 12mm	RC 2128-060 56 / 16mm	
	increased pressure ca. 120g	RC 2124-120 52 / 8mm	RC 2126-120 54 / 12mm	RC 2128-120 56 / 16mm	RC 2130-120 60 / 20mm

Diethrich Non-Traumatic
Cross-Action Bulldog

MICRO BULLDOG



Micro Bulldog
RC 2068
27mm

WELDON (GLOVER) MINI BULLDOG



RC 2069
30mm



RC 2071
28mm



RC 2073
28mm



RC 2076
33mm




RC 2077
31mm




RC 2078
30mm




RC 2079
25mm




RC 2086
25mm



RC 2074
32mm



RC 2080
33mm



RC 2081
25mm



RC 2082
19mm

DEBAKEY NEONATAL ULTRALIGHT CLAMP, 1,2MM



Clamps only available in titanium




RT 2320

DeBakey Neonatal Ultralight Clamp
107mm


A handles straight

















B handles curved




C handles angled



						
C RT 2320 107mm 32mm	B RT 2322 103mm 22mm	B RT 2324 107mm 26mm	A RT 2323 105mm 22mm	A RT 2321 103mm 20mm	B RT 2325 102mm 20mm	B RT 2326 102mm 20mm


						
A RT 2327 111mm 28mm, left	A RT 2328 111mm 28mm, right	A RT 2329 123mm 40mm	B RT 2330 98mm 20mm	B RT 2331 127mm 49mm	A RT 2332 115mm 32mm	

DEBAKEY NEONATAL ULTRALIGHT CLAMP, 1,4MM




DeBakey Ultralight Clamp
jaw tip 1,4mm

Clamps also available in titanium



A	RC 2420	RC 2421	RC 2422	RC 2423	RC 2424	RC 2425	RC 2426
B	RC 2400	RC 2401	RC 2402	RC 2403	RC 2404	RC 2405	RC 2406
C	RC 2440	RC 2441	RC 2442	RC 2443	RC 2444	RC 2445	RC 2446
D	RC 2460	RC 2461	RC 2462	RC 2463	RC 2464	RC 2465	RC 2466



A	RC 2427	RC 2428	RC 2431	RC 2432	RC 2433	RC 2435	RC 2436	RC 2437
B	RC 2407	RC 2408	RC 2411	RC 2412	RC 2413	RC 2415	RC 2416	RC 2417
C	RC 2447	RC 2448	RC 2451	RC 2452	RC 2453	RC 2455	RC 2456	RC 2457
D	RC 2467	RC 2468	RC 2471	RC 2472	RC 2473	RC 2475	RC 2476	RC 2477


COOLEY NEONATAL ULTRALIGHT CLAMP, 1,5MM




A
handles
straight






B
handles
curved



C
handles
angled



Cooley Neonatal Ultralight Clamp
jaw tip 1,5mm



A **RC 6300**
130mm

A **RC 6303**
125mm




A **RC 6306**
115mm

B **RC 6307**
125mm

A **RC 6301**
120mm
B **RC 6304**
117mm

A **RC 6311**
128mm
(6 scale lines)

B **RC 6302**
128mm



A **RC 6299**
120mm

A **RC 6298**
113mm

A **RC 6305**
100mm

B **RC 6308**
117mm

B **RC 6310**
125mm

B **RC 6313**
112mm
(6 scale lines)

C **RC 6296**
120mm


B **RC 6297**
115mm

CASTANEDA NEONATAL CLAMP




RC 6360


Castaneda Neonatal Clamp
handles curved




RC 6360
118mm




RC 6361
116mm




RC 6362
114mm




RC 6316
117mm




RC 6363
105mm




RC 6314
99mm




RC 6315
122mm




RC 6312
118mm




RC 6317
114mm
RC 6365
101mm




RC 6318
116mm




RC 6319
118mm
RC 6366
101mm



RC 6364
115mm



RC 6367
110mm




RC 6368
106mm

DEBAKEY RING HANDLE BULLDOG, 1,2MM



RC 2195

DeBakey Ring Handle Bulldog Clamp
overall length approx. 105-124mm



Clamps also
available in
titanium

RC 2182
134mm
handles straight


RC 2186
130mm
handles straight

RC 2192
100mm
handles straight

RC 2189
122mm
handles curved


RC 2195
129mm
handles curved

MINI-GREGORY CLAMP




Mini-Gregory Clamp
RC 2184
119mm
handles curved

CONTRACTOR BAILEY BABY



THORAX RETRACTOR




Finocchietto Micro Neonatal Thorax Spreader
RS 3050

overall length	overall height	max. spread	blade depth	blade width
59mm	37mm	32mm	10mm	7mm

Art.	overall length	overall height	arm length from top of rack to top	blade depth	blade width	spread
RS 3001	108mm	70mm	55mm	15mm	15mm	75mm
RS 3002	108mm	75mm	60mm	25mm	25mm	75mm
RS 3004	108mm	70mm	60mm	21mm	21mm	75mm

THORAX RETRACTOR


Retractor only available in titanium



Titanium Mini Thorax Spreader
RS 3160SET
 consisting of:
 RS 3160 - spreader frame
 full blade:
 RS 3165 - right
 RS 3166 - left


RS 3160	overall length 76,5mm	overall height 53mm	spread 39mm
----------------	--------------------------	------------------------	----------------

fenestrated blade:
 RS 3167 - right
 RS 3168 - left



full blades: RS 3165 SET

RS 3165 SET			
Art.	blade depth	blade width	type / RIL
RS 3165	15mm	20mm	right
RS 3166	15mm	20mm	left




fenestrated blades: RS 3167 SET

RS 3167 SET			
Art.	blade depth	blade width	type / RIL
RS 3167	15mm	10mm	right
RS 3168	15mm	10mm	left




TITANIUM

THORAX RETRACTOR



overall length	overall height	arm length from top of rack to top	blade depth	blade width	spread
109mm	70mm	56mm	40mm	12mm	65mm




Art.	overall length	overall height	arm length from top of rack to top	blade depth	blade width	spread
RS 3090	81mm	70mm	56mm	15mm	12mm	46mm
RS 3091	109mm	70mm	56mm	15mm	12mm	65mm
RS 3075	127mm	82mm	66mm	20mm	15mm	90mm

SELMAN COOLEY PEDIATRIC CLAMP



Selman Cooley Pediatric Clamp
RC 6170
 120mm
 handles straight




Selman Cooley Pediatric Clamp
RC 6153
 135mm
 handles straight



Selman Cooley Pediatric Clamp
RC 6152
 140mm
 handles straight




RC 6150
 125mm
 handles straight



RC 6151
 120mm
 handles straight

COOLEY PEDIATRIC VASCULAR CLAMP




Cooley Pediatric Vascular Clamp
RC 6520
140mm
handles straight

COOLEY ANASTOMOSIS CLAMP



Cooley Anastomosis Clamp
RC 3125
165mm
handles curved



Cooley Anastomosis Clamp
RC 3845
160mm
handles curved

COOLEY PEDIATRIC CLAMP





RC 6350

Cooley Pediatric Clamp

A
handles
straight




B
handles
curved




A **RC 6350**
145mm
B **RC 6352**
140mm



A **RC 6353**
140mm



A **RC 6351**
140mm




B **RC 6354**
136mm



A **RC 6358**
152mm



A **RC 6571**
142mm




A **RC 6572**
B **RC 6372**
140mm, 20°




A **RC 6573**
B **RC 6373**
138mm



A **RC 6574**
B **RC 6374**
130mm



A **RC 6578**
B **RC 6378**
135mm



A **RC 6579**
131mm


COOLEY PEDIATRIC CLAMP



Clamps also available in titanium



RC 3800
Cooley Patent ductus Clamp




RC 3800
158mm
handles straight

RC 3801
157mm
handles curved




RC 3802
Cooley Coarctation Clamp




RC 3802
170mm
handles straight


RC 3803
168mm
handles curved




RC 3804
Cooley Coarctation Clamp
172mm, handles straight




RC 3805
Cooley Multipurpose Clamp
handles curved



RC 3805
160mm




RC 3806
150mm



RC 3807
130mm


TITANIUM

COOLEY PEDIATRIC CLAMP



RC 5024


Cooley Pediatric Clamp
handles curved




RC 5023
154 / 40mm

RC 5024
159 / 45mm

RC 5025
166 / 52mm



Cooley Peripheral
RC 5026
173mm
handles curved



Cooley Vascular Clamp
RC 5031
160mm
handles straight

COOLEY SPOON-SHAPED CLAMP




Cooley Spoon-Shaped Clamp
RC 3808
180mm
handles straight

COOLEY DERRA CLAMP




RC 3816


Cooley Derra Clamp
handles curved



RC 3815
165mm




RC 3816
165mm




RC 3817
165mm

COOLEY CLAMP




RC 3831
Cooley Derra Clamp
 handles curved

- 1. RC 3830
165mm
- 2. RC 3831
165mm




RC 3835
Cooley Derra Clamp
 handles curved

- 1. RC 3835
175mm
- 2. RC 3836
177mm



RC 3826
Cooley Anastomosis Ligature Clamp
 handles curved

- 1. RC 3826
164mm
- 2. RC 3827
162mm



Cooley Peripheral Vascular External Carotid Clamp
 RC 3839
 185mm, small

KARTCHNER CAROTID CLAMP




Kartchner Carotid Clamp
RC 6401
165 / 56mm
7 scale lines
handles curved

DEBAKEY CLAMP



DeBakey Clamp
RC 6457
145mm
handles angled

COOLEY KITZMILLER PEDIATRIC CLAMP



Cooley Kitzmiller Pediatric Clamp
131mm
handles curved

COOLEY MARTINEZ CLAMP



Cooley Martinez Clamp
RC 6612
180mm, left
handles straight




Cooley Martinez Clamp
RC 6610
180mm, right
handles straight

DEBAKEY PEDIATRIC CLAMP




RC 2196


DeBaKey Pediatric Clamp
145mm
handles curved




RC 2196
14mm
8 scale lines



RC 2197
18mm
9 scale lines



RC 2198
23mm
10 scale lines




RC 2199
25mm
9 scale lines




RC 3135


Derra DeBaKey Clamp
handles curved




DeBaKey Mini Beck
RC 3750
165mm
handles straight



Derra DeBaKey
RC 3135
170mm
handles curved



Derra DeBaKey
RC 3136
172mm
handles curved



Derra DeBaKey
RC 3137
172mm
handles curved

DEBAKEY PEDIATRIC CLAMP




RC 2301


A
handles
straight




B
handles
curved



C
handles
angled



DeBakey Pediatric Clamp




A **RC 2300**
148mm
C **RC 2301**
145mm
B **RC 2310**
145mm

B **RC 2311**
142mm

A **RC 2302**
136mm
C **RC 2314**
130mm
B **RC 2320**
150mm

A **RC 2303**
116mm

A **RC 2306**
137mm
C **RC 2312**
109mm
B **RC 2313**
140mm
B **RC 2315**
145mm



A **RC 2309**
147mm
B **RC 2304**
146mm

A **RC 2307**
147mm

A **RC 2308**
145mm

B **RC 2305**
140mm

DEBAKEY PEDIATRIC CLAMP




Clamps also available in titanium


A handles straight



B handles curved




C handles angled



RC 6454

DeBakey Pediatric Clamp



- A **RC 6440**
163mm
- C **RC 6450**
158mm
- A **RC 6441**
154mm

- C **RC 6451**
158mm


- A **RC 6442**
148mm
- C **RC 6454**
145mm

- A **RC 6443**
130mm

- B **RC 6445**
155mm
- B **RC 6455**
162mm
- B **RC 6453**
155mm
- C **RC 6452**
125mm




- A **RC 6447**
164mm



- B **RC 6444**
165mm



- A **RC 6448**
160mm



- A **RC 6446**
158mm


DEBAKEY PEDIATRIC CLAMP



Model also available in titanium



DeBakey Pediatric Clamp RC 2535
170 / 53mm
handles curved




1. RC 2540
190mm
2. RC 2545
178mm

DeBakey Pediatric Clamp RC 2545
s-form peripheral vascular clamp
handles curved



Dale Clamp RC 2542
handles curved



RC 2541
173mm

RC 2542
178mm

RC 2543
178mm


TITANIUM

DEBAKEY VASCULAR CLAMP



RC 3739


DeBakey Vascular Clamp
handles curved




RC 3735
155mm

RC 3737
148mm

RC 3739
131mm



RC 3731
156mm



RC 3732
158mm

DEBAKEY GLOVER COARCTATION CLAMP



RC 3754

DeBakey Glover Coarctation Clamp



RC 3752
170mm
handles straight



RC 3754
166mm
handles straight


RC 3753
165mm
handles curved

DEBAKEY PERIPHERAL VASCULAR CLAMP




DeBakey Peripheral Vascular Clamp
RC 2550
158mm, handles straight

ERIKSSON-STILLE DEBAKEY CAROTID CLAMP



Eriksson-Stille DeBakey Carotid Clamp
RC 2551
175mm, handles curved

GREGORY PROFUNDA CLAMP



Gregory Profunda Clamp
handles straight


1. RC 2494
145mm, small

2. RC 2495
182mm, medium

3. RC 2496
192mm, large

Gregory Profunda Clamp
RC 2497
115mm
mini-profunda pediatric clamp
handles 50° angled

DEBAKEY BAINBRIDGE CLAMP




Haemostatic Clamp
handles straight

RC 2560
157mm


RC 2565
155mm

RC 2567
180mm

LELAND-JONES CLAMP




Leland-Jones Clamp
handles straight




CASTANEDA PEDIATRIC CLAMP



Castaneda Pediatric Clamp



INSTRUMENTATION SUPPLIES




RC 0800
150mm
2,0mm tip width

RC 0801
150mm
1,5mm tip width

RC 0800

DeBakey Standard Forceps Conical




RC 0809
150mm
1,5mm tip width

RC 0830
150mm
1,0mm tip width

RC 0830

DeBakey Standard Forceps Conical

INSTRUMENTATION SUPPLIES



Baby Mixer Clamp
RA 9494
140mm



Ryder Mini TC-Needle Holder
RN 1810
150mm



Brown Needle Holder
RN 1526
120mm





Brown Needle Holder
RN 1550
130mm



To meet your
specific needs
modifications
are available upon request.





Be sure to check out our online catalog here:
<http://catalog.reuchlen.com>



Expertise
for over
75 years.



August Reuchlen GmbH
Junkersstr. 8
78532 Tuttlingen
Germany
Tel.: +49 (0) 74 61 / 96 52 - 0
Fax: +49 (0) 74 61 / 96 52 - 50
info@reuchlen.com

////////////////////////////////////
www.reuchlen.com